Hygienic Fiber Face Mask Best Practices



Thank you for using a Ruby Lake Hygienic Fiber Face Mask!

This mask has been engineered to provide a top quality filtration layer while staying soft and breathable. As it is a rewashable mask, proper maintenance and cleaning will help it to last as long as possible.

How do I clean my mask?

There are a few different ways you can clean your mask:



HAND WASH

Just like washing your hands, wash the soap in warm, soapy water for 30 seconds



MACHINE WASH

To protect the mask, we recommend you ONLY machine wash the mask in a laundry bag or mesh bag. If you do not have a laundry bag, please do not machine wash mask



SINK OR BASIN

Fill a sink or basin with warm water and a light detergent or soap and wash the mask for at least 30 seconds

Drying Method

Hang dry for 1-2 hours for best results

How many times can I wash my mask?

We engineered the mask to withstand 30-50 wearings total, but with proper care and attention it can last much longer!

Here are things to look out for to know your mask is needing a replacement:

- 1) Loose ear straps
- 2) bagged out, thinning fabric
- 3) loss of tension

30 - 50 WASHINGS

In very dusty environments, the mask may become overloaded with dust or loose particles caught up in the mask and it may cause discomfort. If the dust will not dislodge with a gentle shake or blow of air, then replace the mask with another or wash using the HAND WASH method above to ensure good health and hygiene.